



# WHEELCHAIR SMARTER.

AWARD-WINNING SMART WHEELCHAIR TECHNOLOGY  
RE-IMAGINING MOBILITY AT EVERY TURN.

**LUCI is the first and only smart technology accessory for power wheelchairs, designed to give riders a safer and more inclusive experience.**


LUCI prevents power wheelchair tips, collisions and falls before they happen, thanks to a platform which helps the wheelchair “see” its environment and respond to steering input from the user. The result is safe, confident independence – not autonomy.



Collision  
Avoidance



Drop-off and Tip  
Protection



Cloud-based  
Communication  
and Alerts



Enhanced  
Mobility Platform



Home Assistant  
Compatibility

**LUCI**

LUCI.COM



## LUCI'S 360° VIEW

**LUCI adds hardware and software to riders' new and existing chairs, offering unparalleled stability, security and connectivity.**

Our patented system combines data from sensor fusion, infrared, ultrasonics and radar into a single view of the world, enabling never-before-seen possibilities for power wheelchair riders.

LUCI is currently available for some of the world's top-selling power wheelchairs, including the Permobil F-Series and Permobil M-Series, with standard or alternative access controls. More solutions are coming soon, including compatibility with pediatric power chairs.

SMART WHEELCHAIR TECHNOLOGY FOR  
THE MOST FEARLESS PEOPLE ON EARTH.

## FUELED BY COLLABORATION

We are founded and led by committed members of the community of people living with disabilities and their caregivers. Three years and 17 patents ago, our team set out to make this beautiful, dangerous world more accessible for all people using power wheelchairs. Along the way, we have created a culture of intentional partnership, because technology is another area where inclusion matters.

LUCI and the MyLUCI app are compatible with:

