

# ROVI REQUEST FOR QUOTE

Please complete entire form and send to [orders@luci.com](mailto:orders@luci.com)



## CANADIAN PROVIDER INFORMATION

Provider Company Name:	Billing Address: RECIPIENT NAME:	Shipping Address: RECIPIENT NAME:	<i>Same as billing</i>
Provider Contact Name:	STREET:	STREET:	
Referring Therapist, Clinic, and Contact Info:	CITY:	PROV:	CITY: PROV:
Client Reference / Order Reference Number:	POSTAL CODE:	COUNTRY:	POSTAL CODE: COUNTRY:

## PATIENT CHAIR INFORMATION

Is this for a new wheelchair order or to retrofit an existing wheelchair?  New Chair Order  Retrofit Existing Chair

**PLEASE READ!!** LUCI interface brackets are required for LUCI installation on the ROVI base chairs. If ordering a new chair, select the *LUCI Interface Ready* option on the accessory page of the ROVI form and Motion Concepts will install the brackets at the factory free of charge. If retrofitting to an existing chair, the interface brackets will ship in the box with your LUCI order.

Select LUCI Model: [MSRP: \\$15,051.00 CAD](#)

We do not support swing-away footrests, the extra-large flip-up foot platform, lateral tilt, or high speed motors.

ROVI X3 - LUCI Model 111

ROVI A3 - LUCI Model 112

ROVI A3 Stander - LUCI Model 112

Select Seat-to-Floor Height:

Seat-to-floor height is measured at the front of the seat pan based on a standard seat size. LUCI is not compatible with 17 inch seat to floor height as some seating may interfere with LUCI and/or limit seat travel.

18 inches

19 inches

20 inches

Insurance or Private Pay?

*Pro Tip:* Visit [luci.com/funding/](http://luci.com/funding/) for insurance guidance resources, or connect with our team directly at [education@luci.com](mailto:education@luci.com).

Insurance Carrier

Private Pay

Any custom chair modifications? If so, explain.

Select Drive Method:

This will ensure you receive the correct dashboard mounting bracket:

R-Net LED Remote Joystick

R-Net LED Remote Joystick for Lights

R-Net Color Remote Joystick with Monojacks

R-Net Color Remote Joystick with Monojacks for Lights

R-Net Advanced Joystick

Omni2: ALT Drive Method

Note: We do not support Latch drive or VR2 devices. If you intend to use specialty controls please contact LUCI to confirm compatibility.

Will you be using lower extremity or foot support in addition to the standard footplates (e.g., foot box)?

Yes, please share details in the custom chair modifications box

No

If you are using a configuration not found on this form, please call 615-813-LUCI before ordering.

Send completed form to [orders@luci.com](mailto:orders@luci.com)

ROVI-0426-CAN